Thyroid TSH Rapid Test Cassette

Package Insert
For Self-testing

REF P050502 English

A rapid test for the qualitative detection of human Thyroid Stimulating Hormone (TSH) in human whole blood. For self-testing in vitro diagnostic use only.

[INTENDED USE]
The Thyroid Stimulating Hormone Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of TSH in human whole blood at a cut-off concentration of 5 µIU/mL. This assay provides a preliminary diagnostic test result and can be used for screening for TSH.

[SUMMARY]
Thyroid stimulating hormone (also known as thyrotropin, thyrotropic hormone, TSH, or hTSH for human TSH) is a pituitary hormone that stimulates the thyroid gland to produce thyroxine (T4), and then triiodothyronine (T3) which stimulates the metabolism of almost every tissue in the body.¹ It is a glycoprotein hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland, which regulates the endocrine function of the thyroid.² TSH (with a half-life of about an hour) stimulates the thyroid gland to secrete the hormone thyroxine (T4), which has only a slight effect on metabolism. T4 is converted to triiodothyronine (T3), which is the active hormone that stimulates metabolism. About 80% of this conversion is in the liver and other organs, and 20% in the thyroid itself. Testing of thyroid stimulating hormone³ levels in the blood is considered the best initial test for hypothyroidism. It is important to note the statement from the Subclinical Thyroid Disease Consensus Panel: "There is no single level of TSH at which clinical action is always either indicated or contraindicated. The higher the TSH, the more compelling is the rationale for treatment. It is important to consider the individual clinical context (e.g. pregnancy, lipid profile, APO antibodies)."⁴

The Thyroid Stimulating Hormone Rapid Test Cassette (Whole Blood) is a rapid test that qualitatively detects the presence of TSH in whole blood specimen at the sensitivity of 5 µIU/mL. The Thyroid Stimulating Hormone Rapid Test Cassette (Whole Blood) is a simple test that utilizes a combination of monoclonal antibodies to selectively detect elevated levels of TSH in whole blood.

[PRINCIPLE]
The Thyroid Stimulating Hormone Rapid Test Cassette (Whole Blood) is a qualitative membrane based immunoassay for the detection of Thyroid Stimulating Hormone (TSH) in whole blood. In this test procedure, anti-TSH antibody is immobilized in the test line region and coated particles. After specimen is added to the specimen well of the cassette, it reacts with anti-TSH antibody coated particles in the test. This mixture migrates chromatographically along the length of the test and interacts with the immobilized anti-TSH antibody. Positive specimens react with the specific anti-TSH antibody coated particles to form a colored line at the test line region of the membrane. Absence of this colored line suggests a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

[PRECAUTIONS]
Please read all the information in this package insert before performing the test.

- For self-testing in vitro diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Store in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- This test kit is intended to be used as a preliminary test only and repeatedly abnormal results should be discussed with doctor or medical professional.
- Follow the indicated time strictly.
- Use the test only once. Do not dismantle and touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.



[STORAGE AND STABILITY]
Store as packaged in the sealed pouch at room temperature or refrigerated (2-30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

[MATERIALS PROVIDED]

- Test cassette
- Capillary dropper
- Buffer
- Alcohol pad
- Lancet
- Package insert

[PROCEDURE]

- Wash your hands with soap and rinse with clear warm water.
- Bring the pouch to room temperature before opening it. Open the foil pouch and get out the cassette.
- Carefully pull off the protective released cap of the lancet.
- Use the provided alcohol pad to clean the fingertip of the middle or ring finger as the puncture site.
- Press the lancet, on the side from where the cap was extracted; against the fingertip (Side of ring finger is advised). The tip retracts automatically and safely after use.
- Keeping the hand down massage the end that was pricked to obtain a blood drop.
- Without squeezing the capillary dropper bulb, put it in contact with the blood. The blood migrates into the capillary dropper to the line indicated on the capillary dropper.
- You may massage again your finger to obtain more blood if the line is not reached. Avoid air bubbles.
- Put the blood collected into the sample well of the cassette, by squeezing on the dropper bulb.
- Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and add 2 drops of buffer into the sample well of the cassette.
- Wait for the colored line(s) to appear. Read results at 10 minutes. Do not interpret the result after 20 minutes.

[READING THE RESULTS]

POSITIVE: Two colored lines appear. Both T (Test) and C (Control) line appear.
This result means that the TSH level is higher than the normal (5 µIU/mL) and that you should consult a physician.

NEGATIVE: One colored line appears. Only control line (C) appears.
This result means that the TSH level is not in the range to consider Hypothyroidism.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

[CONTROL PROCEDURE]
A procedural control is included in the test. A colored line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

[LIMITATIONS]

- The Thyroid Stimulating Hormone Rapid Test Cassette (Whole Blood) is for *in vitro* diagnostic use only. The test should be used for the detection of TSH in whole blood specimens only. Neither the quantitative value nor the rate of increase in TSH concentration can be determined by this qualitative test.
- The Thyroid Stimulating Hormone Rapid Test Cassette (Whole Blood) is only for screening the primary hypothyroidism of adult population, not for neonates.

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[BIBLIOGRAPHIE]

- Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
- The American Heritage Dictionary of the English Language, Fourth Edition, Houghton Mifflin Company, 2006. ISBN 0-395-82517-2.
- Sacher R, Richard A, McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company. ISBN 0-8036-0270-7.
- So, M, MacIsaac, R.J., Grossmann M (August 2012), "Hypothyroidism". Australian Family Physician 41 (8): 556-62.
- Surkset al., JAMA 291:228, 2004. Daniel,GH, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, K.J., et al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p. 1666).
- Daniel,GH, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, K.J., et al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p. 1666).

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Lanzette

PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Gangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China

Alkoholtupfer

Ningbo Medsun Medical Co.,Ltd.
No. 55 Jinxi Road, Zhenhai 315221 Ningbo
People's Republic of China

Verteilt von:
Berkeley Health Limited
52 Berkeley Square Mayfair London, W1J 5BT UK
+44 20 8191 0500
support@berkeleyhealth.com www.berkeleyhealth.com

3. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
4. A positive test must be confirmed using a quantitative laboratory TSH assay.
5. False positive results can occur due to heterophile (unusual) antibodies. In certain clinical conditions such as central hypothyroidism, TSH levels may be normal/low, despite hypothyroidism. Medical consultation is recommended to exclude such cases.
6. For Central/Secondary Hypothyroidism, TSH is not a reliable biomarker, which occurs in 1 out of 1,000 Hypothyroidism cases.

【PERFORMANCE CHARACTERISTICS】

An clinical evaluation was conducted comparing the results obtained using the TSH Test to ELISA. The in-house clinical trial included 220 whole blood specimens. The results demonstrated 98.2% specificity and 98.2% sensitivity with an overall accuracy of 98.2%.

Methode	Resultats		ELISA		Total Results
	Positive	Negative	Positive	Negative	
Thyroid Stimulating Hormone Rapid Test Cassette (Whole Blood)	53	163	53	163	56 164
Total Results	54	166	54	166	56 220

Relative Sensitivity: 98.2% (95%CI*: 90.1%-99.9%)
Relative Specificity: 98.2% (95%CI*: 94.8%-99.6%)
Accuracy: 98.2% (95%CI*: 95.4%-99.5%)
*Confidence Interval

【EXTRA INFORMATIONS】

1. How does the TSH test work?
The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the Thyroid Stimulating Hormone Rapid Test Cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3. The capillary dropper provided in the box allows making sure the collected blood volume is correct. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

2. When should the test be used?
The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the Thyroid Stimulating Hormone Rapid Test Cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3. The capillary dropper provided in the box allows making sure the collected blood volume is correct. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

3. Can the result be incorrect?
The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the Thyroid Stimulating Hormone Rapid Test Cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3. The capillary dropper provided in the box allows making sure the collected blood volume is correct. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

4. How to interpret the test if the color and the intensity of the lines are different?
The color and intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. The test should be considered as positive whatever the color intensity of the test line is.

5. If I read the result after 20 minutes, will the result be reliable?
No. The result should be read at 10 minutes after adding the buffer. The result is unreliable after 20 minutes.

6. What do I have to do if the result is positive?
The result is positive, it means that the TSH level in blood is higher than the normal (5 µIU/mL) and that you should consult a physician to show the test result. Then, the physician will decide whether additional analysis should be performed.

7. What do I have to do if the result is negative?
If the result is negative, it means that the TSH level is below 5 µIU/mL and is within the normal range. A case of Hypothyroidism, though rare, but cannot be excluded based on such test results. However, if the symptoms persist, it is recommended to consult a physician.

【BIBLIOGRAPHY】

- Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
- The American Heritage Dictionary of the English Language, Fourth Edition, Houghton Mifflin Company, 2006. ISBN 0-395-82517-2.
- Sacher R, Richard A, McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company. ISBN 0-8036-0270-7.
- So, M, MacIsaac, R.J., Grossmann M (August 2012), "Hypothyroidism". Australian Family Physician 41 (8): 556-62.
- Surkset al., JAMA 291:228, 2004. Daniel,GH, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, K.J., et al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p. 1666).

Index of symbols		Tests per kit		Authorized representative in EU	
	Consult instructions for use		Use by		Representant autorisé dans l'UE
	For <i>in vitro</i> diagnostic use only		Lot number		Ne pas réutiliser
	Store between 2-30 °C		Manufacturer		Réf. catalogue
	Do not use if package is damaged				

Hangzhou AllTest Biotech Co.,Ltd.
#550, Yinhai Street
Hangzhou Economic & Technological Development Area
Hangzhou, 310118 P.R. China
Web: www.alltests.com.cn Email: info@alltests.com.cn

Lancet

PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Gangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China

Alcohol Pad

Ningbo Medsun Medical Co.,Ltd.
No. 55 Jinxi Road, Zhenhai 315221 Ningbo
People's Republic of China

Shanghai International Holding Corp. GmbH (Europe)
Eiffelstraße 80
20537 Hamburg, Germany

Distributed by:
Berkeley Health Limited
52 Berkeley Square Mayfair London, W1J 5BT UK
+44 20 8191 0500
support@berkeleyhealth.com www.berkeleyhealth.com

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[BIBLIOGRAFÍA]

- Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
- The American Heritage Dictionary of the English Language, Fourth Edition. Houghton Mifflin Company.2006.ISBN 0-395-82517-2.
- Sacher R, Richard A, McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company.ISBN 0-8036-0270-7.
- So, M, MacIsaac, RJ, Grossmann M (August 2012). "Hypothyroidism". *Australian Family Physician* **41** (8): 656-62.
- Surkesl.at. JAMA 291:228, 2004. Daniel.G.H, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, KJ, et al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p.1666).

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[BIBLIOGRAFIA]

- Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
- The American Heritage Dictionary of the English Language, Fourth Edition. Houghton Mifflin Company. 2006. ISBN 0-395-82517-2.
- Sacher R, Richard A, McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company.ISBN 0-8036-0270-7.
- So, M, MacIsaac, RJ, Grossmann M (August 2012). "Hypothyroidism". *Australian Family Physician* **41** (8): 556-62.
- Surkesl.at. JAMA 291:228, 2004. Daniel.G.H, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, KJ, et al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p. 1666).

Hangzhou AllTest Biotech Co.,Ltd. #550, Yinhai Street Hangzhou Economic & Technological Development Area Hangzhou, 310018 P.R., China Web: www.alltests.com.cn Email: info@alltests.com.cn			

PROMISEMED HANGZHOU MEDITECH CO., LTD. No. 1388 Cangxing Street, Cangqian Community, Yuhang District, Hangzhou City, 311121 Zhejiang, China			

Ningbo Medsun Medical Co., Ltd. No. 55 Jinx Road, Zhenhai 315221 Ningbo People's Republic of China			

Distribuito da: Berkeley Health Limited 52 Berkeley Square Mayfair London, W1J 5BT UK +44 20 8191 0500 support@berkeleyhealth.com www.berkeleyhealth.com			

Hangzhou AllTest Biotech Co.,Ltd. #550, Yinhai Street Hangzhou Economic & Technological Development Area Hangzhou, 310018 P.R., China Web: www.alltests.com.cn Email: info@alltests.com.cn			

PROMISEMED HANGZHOU MEDITECH CO., LTD. No. 1388 Cangxing Street, Cangqian Community, Yuhang District, Hangzhou City, 311121 Zhejiang, China			

Ningbo Medsun Medical Co., Ltd. No. 55 Jinx Road, Zhenhai 315221 Ningbo People's Republic of China			

--	--	--	--