**Urinary Tract Infections Rapid Test  
Package Insert  
For Self-testing**

REF P050515 English

For rapid qualitative detection of Leukocytes, Blood, Nitrite and Protein in human urine.

For self-testing in vitro diagnostic use only.

**INTENDED USE**

The Urinary Tract Infections Test (Urine) is firm plastic strips onto which several separate reagent areas are affixed. The test is for the qualitative detection of the following analytes in urine: Leukocytes, Blood, Nitrite and Protein. The Urinary Tract Infections Test (Urine) is for single use in self-testing.

**SUMMARY**

A urinary infection represents the most common disease of the urinary tract which includes the urethra, the bladder, the ureter and the kidneys. Men, women and children are likely to experience a urinary infection. It's mostly women who suffer from urinary infections, since the short urethra favours the penetration of germs. However, elderly males are also affected if they have an enlarged prostate which obstructs the urine flow.

In healthy people, urine is sterile (i.e. it doesn't contain any micro-organisms). One of the best ways to keep your urinary tract sterile is to empty your bladder completely at regular intervals. Generally, an infection starts in the urethra and may then spread into the upper urinary tract as far as the kidneys.

The symptoms vary considerably: burning when emptying the bladder, or a strong urge to urinate. The urine may also be cloudy or have a strong odour.

**PRINCIPLES OF THE EXAMINATION METHOD**

Leukocytes: This test reveals the presence of granulocyte esterases. The esterases cleave a derivatized pyrazole amino acid ester to liberate derivatized hydroxyl pyrazole. This pyrazole then reacts with a diazonium salt to produce a beige-pink to purple color.

Blood: This test is based on the peroxidase-like activity of hemoglobin which catalyzes the reaction of disopropylbenzene dihydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting color ranges from orange to green to dark blue.

Nitrite test depends upon the conversion of nitrate to nitrite by the action of Gram-negative bacteria in the urine. In an acidic medium, nitrite in the urine reacts with p-arsanilic acid to form a diazo-like compound. This diazonium compound in turn couples with 1-(N-(1-naphthyl)ethyl)enediamine to produce a pink color.

Protein: This reaction is based on the phenomenon known as the "protein error" of pH indicators (Tetrabromophenol Blue). The anion produced by pH indicators under given conditions combines with the cation produced by protein, then pH indicators colors from yellow to green-blue for positive results.

**PRECAUTIONS****Read the instructions carefully before performing the test.**

• For self-testing in vitro diagnostic use only.

• Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.

• A clean container - uncontaminated by cleaning fluids - to collect urine.

• Keep out of the reach of children.

• Do not use after the expiry date or if the pouch is damaged.

• Follow the indicated time strictly.

• Use the test only once. Do not dismantle and touch the reagent areas of the test strip.

• For self-testing only.

• The used test should be discarded according to local regulations.

• In case of difficulties in color identification (such as Daltonism), ask for help in test reading.

**STORAGE AND STABILITY**Store as packaged at room temperature or refrigerated (2-30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE**. Do not use beyond the expiration date.**MATERIALS PROVIDED**

• Test strip • Plastic cup • Color chart • Package insert

• Specimen container

**MATERIALS REQUIRED BUT NOT PROVIDED**

• Timer or a watch with a second hand

**PROCEDURE**

ATTENTION: It's recommended to take a sample of urine for the test in the early morning since it's the most concentrated. The urine used for the test should not come into contact with water from the toilet or any disinfectant or cleaning substances.

For women only: The test should not be performed during or for three days after your menstrual period. The urine sample should not be contaminated with vaginal fluids since this may produce a misleading result.

Do not make any important medical decision without first referring to your doctor.

**COLLECT URINE:**


Collect part of the urine in the supplied plastic cup or using a clean cup without any residual detergents. Make sure to fill up the cup with urine.

**PERFORMING THE TEST:**1) Open the foil pouch and take out the test strip. **Do not touch the test fields**. Once opened the pouch, it is recommended to perform the test **immediately**.

2) Dip the test strip in the urine sample.

ATTENTION: Press the strip and make sure that all four test fields are immersed for about **1-2 seconds**.

3) Then remove the test strip and wipe off any surplus urine against the rim of the container or with an absorbent material (e.g. a paper towel) to avoid mixing chemicals from adjacent reagent areas.

4) Wait for **2 min.** (do not read results after 3 minutes) Read the result separately for each parameter, compare color with color chart provided.**READING THE RESULTS**

Read the result separately for each parameter; compare color with color chart provided.

Color changes on the edges of the test pads or color changes after more than 3 minutes have to be ignored.

**NEGATIVE**  
The test field for LEUKOCYTES stayed whitish.  
The Test field for BLOOD stayed mustard yellow.  
The Test field for NITRITE stayed white.  
The Test field for PROTEIN stayed yellowish.

**POSITIVE RESULT FOR LEUKOCYTES**  
If the color of the test field has changed to **purple**, then leukocytes have been found in your urine.

**POSITIVE RESULT FOR BLOOD**  
If the color of the test field has changed to **green** (or some green spots appear on the background), then blood has been found in your urine.

**POSITIVE RESULT FOR NITRITE**  
If the color of the test field has changed to **pink**, then nitrites have been found in your urine.

**POSITIVE RESULT FOR PROTEIN**  
If the color of the test field has changed to **green**, then proteins have been found in your urine.

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