


berkeleyhealth DRUG

Self-Test for the detection of 6 drugs in urine: Morphine, Amphetamine, Cocaine, Marijuana, Methamphetamine and Methadone



CE
0483




DETECTED SUBSTANCES

TEST	Cut-off (ng/mL)	Time permanence in urine after exposure
Amphetamine (AMP1000)	1000	4 to 24 hours
Cocaine (COC 300)	300	24 to 48 hours
Marijuana	50	3 to 10 days
Methadone (MTD 300)	300	3 to 5 days
Methamphetamine (MET 500)	500	9 to 24 hours
Morphine (MOP 300)	300	3 to 5 days

TECH SPECS

	CUT-OFF	SENSITIVITY	SPECIFICITY	OVERALL ACCURACY
AMP	1000 ng/mL	99.9%	99.9%	99.9%
MET	500 ng/mL	99.9%	99.9%	99.9%
COC	300 ng/mL	99.9%	99.9%	99.9%
MOP	300 ng/mL	99.9%	99.9%	99.9%
MTD	300 ng/mL	99.9%	99.9%	99.9%
THC	50 ng/mL	99.9%	99.9%	99.9%



WHO ARE THE INTENDED USERS

DRUG TEST can be used by everyone who wants to check the presence of the above mentioned substances in urine.

TEST PRINCIPLE

DRUG TEST is an immunochromatographic assay that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect levels of specific drugs in urine.

CONTENT:

1 sealed aluminium pouch containing: 1 test panel and 1 desiccant bag; 1 instructions for use leaflet.

HOW TO USE IT

- 1) Open the cap, with the arrow pointing toward the urine specimen, dip the test panel vertically in the urine specimen for at least 10 to 15 seconds.
- 2) Replace the cap and leave the test panel on a non-absorbent flat surface.
- 3) Wait 5 minutes before reading the result.