Distributed by:
berkeleyhealth
52 Berkeley Square
Mayfair London W1J 5BT UK
+44 20 8191 0500
support@berkeleyhealth.com
www.berkeleyhealth.com

EC REP Qarad EC-REP BV
Pas 257
2440 Geel - Belgium

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
+41 91 605 1030
support@primallabs.ch
primallabs.ch

REF 890080

890080_IFU_72_5.0 08/2022

Registered trademarks.

berkeleyhealth PROSTATE PSA

PROSTATE PREVENTION

- Self-test for the determination of the Prostate-Specific Antigen (PSA) in whole blood samples
- Selbsttest zum Nachweis des Prostate-spezifischen Antigens (PSA) in Vollblutproben
- Auto-test pour la détection de l'Antigène Prostatique Spécifique (PSA) dans des échantillons de sang entier
- اختبار ذاتي لتحديد اجمالي مستضد البروستاتا النوعي في عينات الدم الكامل

INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
NOTICE D'UTILISATION
نشرة تعليمات الاستخدام

PROSTATE AND PSA

The Prostate Specific Antigen (PSA) is a protein produced by the epithelial cells of the prostate gland. Its function is to fluidify the ejaculate and therefore increase sperm mobility. PSA also plays a role in dissolving cervical mucus, allowing sperm to enter the uterus. Small quantities of PSA are normally found in the blood of men with a healthy prostate. Prostate pathologies (such as prostatitis, benign prostatic hyperplasia, cancer) lead to an increase in the blood levels of PSA. PSA levels also increase naturally with age, after prostate surgery or after specific diagnostic tests (rectal examination, prostate biopsy, prostate ultrasound scan, etc.).

WHO IS THE PROSTATE PSA TEST FOR?


- For men who do not have symptoms, are aged over 45, who wish to check their prostate;
- For men with urinary tract symptoms (difficulty urinating, pain, increased frequency).

REFERENCES

- Catalona WJ, Southwick PC, Slawin KM, Partin AW, Brawer M.K., Flanigan R.C., Patel A., Richie J.P., Walsh P.C., Scardino P.T., Lange P.H., Gasior G.H., Loveland K.G., Bray K.R.: Comparison of percent free PSA, PSA density, and age-specific PSA cutoffs for prostate cancer detection and staging. Urology. 2000 Aug 1;56(2):255-60.
- Fritsche H.A., Babaian R.J.: Analytical performance goals for measuring prostate-specific antigen. Clin Chem. 1993 Jul;39(7):1525-29.
- Catalona WJ. "Prostate Cancer Screening" Med Clin North Am. 2018 Mar;102(2):199-214.

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- * Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
- * Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.
- * The desiccant packet must not be used. Dispose of it with household waste without opening it.

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







TEST PROCEDURE

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag.
3) Carefully rotate the protective cap of the sterile lancet 360° without pulling it.
4) Carefully massage the finger chosen for the puncture...
5) Holding the hand down, massage the finger until a large drop of blood forms.
6) Take the pipette without pressing the bulb.
7) Place the blood collected with the pipette into the well indicated on the cassette (S) by pressing the pipette bulb.
8) Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on).
9) Wait 5 minutes and read the results referring to the next section for result interpretation.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES
The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line. This result means that the value of Prostate-Specific Antigen (PSA) in the blood is higher than normal, so you must consult a physician.

NEGATIVE RESULTS

A coloured band appears only under the C (control) sign. This means that the value of Prostate-Specific Antigen (PSA) in the blood is normal. A negative result does not exclude the possibility of prostatic pathologies.

NON VALID RESULT

No bands appear or there is a line only under the T (test) sign and not under the C (control) sign. In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test with a new blood sample.

FAQ - QUESTIONS AND ANSWERS

HOW DOES THE PROSTATE PSA TEST WORK? The Prostate-Specific Antigen is a protein secreted normally by the prostate and circulating in our body in small quantities. Values above normal may be a sign of the presence of pathological conditions in the prostate, which may be of infectious or inflammatory origin...
WHEN CAN THE TEST BE USED? The PROSTATE PSA Test can be performed at any time of the day, but in some cases it is advisable to wait for a specific period of time in order to avoid a false positive result.
CAN THE RESULT BE INCORRECT? The result is correct if the instructions are followed carefully. However, the result may not be correct if: the device comes into contact with other liquids before being used, if part of

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